

APP669-711

## Chemical Properties

CAS No. :

Formula: C210H325N57O61S2

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	APP669-711 is a peptide fragment comprising amino acids 669 to 711 of the amyloid precursor protein (APP). This peptide serves as a biomarker in Alzheimer's disease (AD) research and is useful for diagnosing amyloid deposition in the brain.
Targets(IC50)	Beta Amyloid

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481