

BCOR(498-514), biotinylated

Chemical Properties

CAS No. :

Formula: C98H155N25O29S

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	BCOR(498-514), biotinylated, is the smallest BCL6 binding domain with a KD of 1.32 μ M. It can block BCL6-mediated transcriptional repression and induce cell death in lymphoma cells.
Targets(IC50)	Bcl-6

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481