

Maleimide-PEG8-Val-Ala-PAB

Chemical Properties

CAS No. :

Formula: C41H65N5O15

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Maleimide-PEG8-Val-Ala-PAB is an ADC linker composed of a cleavable Val-Ala-PAB linker and the probe molecule Maleimide. It is used to connect the monoclonal antibody Cetuximab with the CDK inhibitor SNS-032 for the synthesis of antibody-drug conjugates (ADC).
Targets(IC50)	ADC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481