

Highly purified Type V collagen, from human placenta

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Type V collagen, highly purified from human placenta (Human Placenta Type V collagen, immunization grade), is a form of collagen protein that induces an immune response in animals, leading to the production of specific antibodies against this collagen. It also serves as a substrate for matrix metalloproteinases (MMP).
Targets(IC50)	MMP

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481