

Alkyne-C6-OMs

Chemical Properties

CAS No. : 228403-45-0


Formula: C₉H₁₆O₃S

Molecular Weight: 204.29

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Alkyne-C6-OMs is a PROTAC linker utilized in the synthesis of PROTAC compounds.
Targets(IC50)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.895 mL	24.475 mL	48.950 mL
5 mM	0.979 mL	4.895 mL	9.790 mL
10 mM	0.4895 mL	2.4475 mL	4.895 mL
50 mM	0.0979 mL	0.4895 mL	0.979 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481