

Highly purified Type V collagen, from human amnion

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Human Amnion Type V collagen (immunization grade) is a highly purified form of collagen derived from human amnion. It has the ability to stimulate the immune system of animals to produce antibodies specific to this collagen. Additionally, this collagen serves as a hydrolytic substrate for MMP (matrix metalloproteinases).
Targets(IC50)	MMP

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481