

## L-Lactate Dehydrogenase, bovine muscle

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	L-Lactate Dehydrogenase, bovine muscle, is a redox enzyme that catalyzes the reversible conversion of lactate to pyruvate while simultaneously reducing NAD <sup>+</sup> to NADH.
Targets(IC50)	Endogenous Metabolite

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481