

Highly purified Type III collagen, from human placenta

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Human Type III collagen, derived from human placenta (immunization grade), is a highly purified compound capable of eliciting an immune response in animals to produce specific antibodies against this collagen protein. Additionally, this collagen serves as a substrate for the hydrolysis by matrix metalloproteinases (MMPs).
Targets(IC50)	MMP

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481