

Human Kallikrein

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	HumanKallikrein is a serine protease found in both plasma and tissues. It has the potential to be studied for its roles in regulating blood pressure, activating the complement system, and mediating and sustaining inflammatory responses.
Targets(IC50)	Serine Protease

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481