

Sulfo-Cyanine7 NHS ester TEA

Chemical Properties

CAS No. :

Formula: C47H62N4O10S2

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sulfo-Cyanine7 NHS ester (TEA) is an amine-reactive succinimidyl ester. This reagent allows for the straightforward preparation of sulfo-Cyanine7-labeled biomolecules, such as proteins. The dye-labeled molecules can then be utilized in a variety of experiments related to research and active molecule design.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481