

[D-His26]-Neuropeptide Y, human, rat

Chemical Properties

CAS No. :

Formula: C180H276N54O55S

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	[D-His26]-Neuropeptide Y, human, rat is a synthetic variant of neuropeptide Y (NPY). This compound acts as a Y1R agonist and may help prevent the development of anxiety, social disorders, and depressive symptoms. It has potential as an early intervention treatment for post-traumatic stress responses.
Targets(IC50)	Neuropeptide Y Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481