

Highly purified Type I collagen, from human placenta

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Human Type I collagen, derived from highly purified human placenta and classified as immunization grade, stimulates the animal immune system to produce specific antibodies against this collagen. Additionally, collagen serves as a hydrolytic substrate for MMP (Matrix Metalloproteinases).
Targets(IC50)	MMP

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481