

## 5-Lipoxygenase blocking peptide

### Chemical Properties

CAS No. :

Formula: C122H188N42O33S

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	5-Lipoxygenase blocking peptide (Human/rat 5-LO 130-149) is a specific sequence fragment of 5-lipoxygenase (5-LOX) used for generating antibodies against 5-LOX.
Targets(IC50)	Lipoxygenase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481