

## Benzyl-PEG5-THP

## Chemical Properties

CAS No. : 230620-74-3


Formula: C<sub>22</sub>H<sub>36</sub>O<sub>7</sub>

Molecular Weight: 412.52

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Benzyl-PEG5-THP is a PROTAC linker, classified under the PEG category. It is utilized in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4241 mL	12.1206 mL	24.2412 mL
5 mM	0.4848 mL	2.4241 mL	4.8482 mL
10 mM	0.2424 mL	1.2121 mL	2.4241 mL
50 mM	0.0485 mL	0.2424 mL	0.4848 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481