

## Benzyl-PEG8-MS

## Chemical Properties

CAS No. : 1330059-91-0


Formula: C<sub>24</sub>H<sub>42</sub>O<sub>11</sub>S

Molecular Weight: 538.65

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Benzyl-PEG8-MS is a PEG-based PROTAC linker used in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8565 mL	9.2825 mL	18.5649 mL
5 mM	0.3713 mL	1.8565 mL	3.713 mL
10 mM	0.1856 mL	0.9282 mL	1.8565 mL
50 mM	0.0371 mL	0.1856 mL	0.3713 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481