

Bromo-PEG3-THP

Chemical Properties

CAS No. : 220568-33-2


Formula: C₁₁H₂₁BrO₄

Molecular Weight: 297.19

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Bromo-PEG3-THP is a PEG-based PROTAC linker, used in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3649 mL	16.8243 mL	33.6485 mL
5 mM	0.673 mL	3.3649 mL	6.7297 mL
10 mM	0.3365 mL	1.6824 mL	3.3649 mL
50 mM	0.0673 mL	0.3365 mL	0.673 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481