


t-Boc-N-amido-PEG2-C6-Cl

Chemical Properties

CAS No. :	920264-35-3
Formula:	C ₁₅ H ₃₀ ClNO ₄
Molecular Weight:	323.86
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	t-Boc-N-amido-PEG2-C6-Cl is a PROTAC linker of the PEG class, utilized in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0878 mL	15.4388 mL	30.8775 mL
5 mM	0.6176 mL	3.0878 mL	6.1755 mL
10 mM	0.3088 mL	1.5439 mL	3.0878 mL
50 mM	0.0618 mL	0.3088 mL	0.6176 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481