


Benzoylarginine leucinamide

Chemical Properties

CAS No. :	66127-57-9
Formula:	C ₁₉ H ₃₀ N ₆ O ₃
Molecular Weight:	390.48
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Benzoylarginine leucinamide is a dipeptide.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.561 mL	12.8048 mL	25.6095 mL
5 mM	0.5122 mL	2.561 mL	5.1219 mL
10 mM	0.2561 mL	1.2805 mL	2.561 mL
50 mM	0.0512 mL	0.2561 mL	0.5122 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Blanco RM, Alvaro G, Guisán JM. Enzyme reaction engineering: design of peptide synthesis by stabilized trypsin. Enzyme Microb Technol. 1991 Jul;13(7):573-83. PubMed PMID: 1367640.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481