

## Free Aliphatic Acid Assay Kit (Spectrophotometry)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years  
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Free aliphatic acids (Free Fatty Acids, FFA) are both products of lipid hydrolysis and substrates for lipid synthesis. The concentration of FFA is associated with lipid metabolism, glucose metabolism, and endocrine function, and can also reflect quality changes during food storage.
-------------	--

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481