

## Triglyceride Content Assay Kit (Spectrophotometry)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years  
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Triglycerides (TG) are lipid molecules formed from long-chain fatty acids and glycerol. They are not only major components of cell membranes but also serve as important respiratory substrates.
-------------	--

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481