

Triglyceride Content Assay Kit (Microanalysis)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at low temperature,Keep away from direct sunlight

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Triglycerides (TG) are lipid molecules formed from long-chain fatty acids and glycerol. They are not only major components of cell membranes but also serve as important respiratory substrates.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481