

Compound N023-0008

## Chemical Properties


CAS No. :

Formula: C<sub>29</sub>H<sub>43</sub>N<sub>05</sub>

Molecular Weight: 485.665

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound N023-0008 is a useful organic compound for research related to life sciences and the catalog number is T131314.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.059 mL	10.2951 mL	20.5901 mL
5 mM	0.4118 mL	2.059 mL	4.118 mL
10 mM	0.2059 mL	1.0295 mL	2.059 mL
50 mM	0.0412 mL	0.2059 mL	0.4118 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481