


Pulchinenoside J

Chemical Properties

CAS No. :	1476801-10-1
Formula:	C ₆₅ H ₁₀₆ O ₃₀
Molecular Weight:	1367.52
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pulchinenoside J is a natural product suitable for biochemical experiments and drug synthesis research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7313 mL	3.6563 mL	7.3125 mL
5 mM	0.1463 mL	0.7313 mL	1.4625 mL
10 mM	0.0731 mL	0.3656 mL	0.7313 mL
50 mM	0.0146 mL	0.0731 mL	0.1463 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481