

Anti-Butyryl-Histone H2B (Lys12) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Purification:	Antigen affinity purification

Applications

Verified Activity:	Blocking buffer: 5% NFD/MTBST Primary ab dilution: 1: 2000 Primary ab incubation condition: 2 hours at room temperature Secondary ab: Goat Anti-Rabbit IgG H&L (HRP) Lysate: (-) HeLa, (+) HeLa+Sodium butyrate (30 mM, 4 h) Protein loading quantity: 20 µg Exposure time: 60 s Predicted MW: 14 kDa Observed MW: 14 kDa
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Research Background

Histone post-translational modifications (PTMs) are key mechanisms in epigenetic regulation of chromatin structure, often referred to as the "histone code." PTMs on histones include acetylation, methylation, phosphorylation, and several novel acylation modifications discovered in recent years. These histone modifications directly influence chromatin structure and the binding of transcription factors or other epigenetic regulators, thereby altering genomic stability and gene transcription. Histone methylation typically occurs on lysine and arginine residues of core histones. It can either promote or inhibit gene transcription, depending on whether methylation occurs on lysine or arginine residues and the number of methyl groups present (lysines can be mono-, di-, or trimethylated, while arginines can undergo mono-, symmetric, or asymmetric dimethylation). Lysine methylation in histones commonly occurs at lysine residues 4, 9, 27, 36, and 79 of histone H3, and at lysine 20 of histone H4. Arginine methylation typically occurs at arginine residues 2, 8, 17, and 26 of histone H3, and at arginine 3 of histone H4. Histone methyltransferases (HMTs) and histone demethylases (HDMs) are the primary regulators of methylation modifications.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481