

## Anti-CDX1 Antibody (1U448)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 28 kDa. Actual: 28 kDa.
Clone:	1U448
Purification:	Protein A purified

### Applications

Verified Activity:	Western blot analysis of Caco-2 cell lysate. Using CDX1 (TMAB-04192) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB,ICC/IF
Recommended	WB: 1:1000-2000; ICC/IF: 1:50-200

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human CDX1
Antigen Species:	Human
Gene ID:	1044
Uniprot ID:	P47902

### Research Background

Could play a role in the terminal differentiation of the intestine.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481