

Anti-Coilin Antibody (9W479)

Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 63 kDa. Actual: 76 kDa.
Clone:	9W479
Purification:	Protein A purified

Applications

Verified Activity:	Western blot analysis of Jurkat cell lysate. Using Coilin (TMAB-04513) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB,FCM
Recommended	WB: 1:1000-2000; FCM: 1:50-100

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Coilin
Antigen Species:	Human
Gene ID:	8161
Uniprot ID:	P38432

Research Background

Component of nuclear coiled bodies, also known as Cajal bodies or CBs, which are involved in the modification and assembly of nucleoplasmic snRNPs.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481