

Anti-DAZAP1 Antibody (9R964)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Molecular Weight:	Theoretical: 43 kDa. Actual: 43 kDa.
Clone:	9R964
Purification:	Protein A purified

Applications

Verified Activity:	<p>1. Western blot analysis of Jurkat cell lysate. Using DAZAP1 (TMAB-04990) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.</p> <p>2. 4% Paraformaldehyde-fixed Jurkat (H) cell; Triton X-100 at r. T. for 20 min; Antibody incubation with (DAZAP1) monoclonal Antibody, unconjugated (TMAB-04990) 1:100, 90 min at 37°C; followed by conjugated Goat Anti-Rabbit IgG antibody (green) at 37°C for 90 min, DAPI (blue) was used to stain the cell nuclei. PBS instead of the primary antibody was used as the blank control.</p> <p>3. The Jurkat (H) cells were fixed with 4% PFA (10 min at r. T.) and then permeabilized with 90% ice-cold methanol for 20 min at -20°C, the cells then were incubated in 5% BSA to block non-specific protein-protein interactions (30 min at r. T.), followed by secondary antibody incubation for 40 min at room temperature. Primary Antibody (green): Rabbit Anti-DAZAP1 antibody (TMAB-04990, 1:100); Isotype Control (orange): Rabbit IgG. Blank control (black): PBS. Acquisition of 20,000 events was performed.</p>
Application:	WB,FCM,ICC/IF
Recommended	WB: 1:1000-2000; FCM: 1:50-100; ICC/IF: 1:50-200

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human DAZ associated protein 1
Antigen Species:	Human
Gene ID:	26528
Uniprot ID:	Q96EP5

Research Background

RNA-binding protein, which may be required during spermatogenesis.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481