

Anti-DcR2 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 36 kDa. Actual: 50 kDa.
Purification:	Protein A purified

Applications

1. Blank control: Jurkat cells (blue). Primary Antibody: Rabbit Anti-DcR2 antibody (TMAB-05012), Dilution: 1 μ g in 100 μ L 1X PBS containing 0.5% BSA; Isotype Control Antibody: Rabbit IgG (orange), used under the same conditions); Secondary Antibody: Goat anti-rabbit IgG-PE (white blue), Dilution: 1: 200 in 1 X PBS containing 0.5% BSA.

Protocol

The cells were fixed with 2% paraformaldehyde (10 min). Primary antibody (TMAB-05012, 1 μ g /1x10⁶ cells) were incubated for 30 min on the ice, followed by 1 X PBS containing 0.5% BSA+10% goat serum (15 min) to block non-specific protein-protein interactions. Then the Goat Anti-rabbit IgG/PE antibody was added into the blocking buffer mentioned above to react with the primary antibody at 1/200 dilution for 30 min on ice. Acquisition of 20,000 events was

Verified Activity:

performed.

2. Sample:

Lane 1: Human MKN45 cell lysates

Lane 2: Human SW480 cell lysates

Lane 3: Human HUVEC cell lysates

Lane 4: Human DU145 cell lysates

Lane 5: Human HepG2 cell lysates

Primary: Anti-DcR2 (TMAB-05012) at 1/1000 dilution

Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution

Predicted band size: 36 kDa

Observed band size: 50 kDa

Application:

WB,FCM

Recommended

WB: 1:500-2000; FCM: 1 μ g/Test

Properties

Stability & Storage:

Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping:

Shipping with blue ice.

Antigen Details

Immunogen: KLH conjugated synthetic peptide: human DcR2
Antigen Species: Human
Gene ID: 8793
Uniprot ID: Q9UBN6

Research Background

Apoptosis is induced by certain cytokines including TNF and Fas ligand in the TNF family through their death domain containing receptors. TRAIL/Apo2L, a member of the TNF family, induces apoptosis of a variety of tumor cell lines. DR4 and DR5 are functional receptors for TRAIL, and DcR1/TRID is a decoy receptor. Another member of the TRAIL receptor family was identified and designated DcR2. The DcR2 receptor is 386 amino acids in length and has an extracellular TRAIL binding domain, but lacks intracellular death domain and does not induce apoptosis. Although this receptor binds to the cytotoxic ligand TRAIL, it contains a truncated death domain and functions as an inhibitory receptor. When overexpressed, the DcR2 receptor can protect cells against TRAIL mediated cytotoxicity. Like DR4 and DR5, DcR2 transcript is widely expressed in a variety of normal human tissues but DcR2 is absent in most tumors. Ultraviolet radiation has been shown to upregulate DcR2 expression on human keratinocytes. Over expression of DcR2 attenuated TRAIL induced apoptosis.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481