

## Anti-FAM120B Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Mouse,Rat,Cow,Horse)
Molecular Weight:	Theoretical: 103 kDa. Actual: 80 kDa.
Purification:	Protein A purified

### Applications

Sample:	HL60 (Human) Cell Lysate at 30 µg U2os (Human) Cell Lysate at 30 µg
Verified Activity:	Primary: Anti-FAM120B (TMAB-05820) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 80/93 kD Observed band size: 80 kD
Application:	WB
Recommended	WB: 1:500-2000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human FAM120B
Antigen Species:	Human
Gene ID:	84498
Uniprot ID:	Q96EK7

### Research Background

Functions as a transactivator of PPARG and ESR1. Functions in adipogenesis through PPARG activation.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481