

## Anti-FAM220A Polyclonal Antibody

### Product Details

|                   |                      |
|-------------------|----------------------|
| Ig Type:          | IgG                  |
| Reactivity:       | Human                |
| Molecular Weight: | Theoretical: 28 kDa. |
| Purification:     | Protein A purified   |

### Applications

|                    |   |
|--------------------|---|
| Verified Activity: | Paraformaldehyde-fixed, paraffin embedded (human gastric carcinoma); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30 min; Antibody incubation with (FAM220A) Polyclonal Antibody, Unconjugated (TMAB-05859) at 1:400 overnight at 4°C, followed by a conjugated secondary for 20 minutes and DAB staining. |
| Application:       | IHC-P,IHC-Fr,IF   |
| Recommended        | IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500  |

### Properties

|                      |   |
|----------------------|---|
| Stability & Storage: | Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping:            | Shipping with blue ice.   |

### Antigen Details

|                  |   |
|------------------|---|
| Immunogen:       | KLH conjugated synthetic peptide: human FAM220A |
| Antigen Species: | Human   |
| Gene ID:         | 84792   |
| Uniprot ID:      | Q7Z4H9  |

### Research Background

May negatively regulate STAT3. Interacts with STAT3.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481