

## Anti-FOXD4L2 Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 46 kDa. Actual: 46 kDa.
Purification:	Protein A purified

### Applications

Verified Activity:	Sample: JurkaT (Human) Cell Lysate at 30 µg Primary: Anti-FOXD4L2 (TMAB-06119) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 46 kD Observed band size: 46 kD
Application:	WB
Recommended	WB: 1:500-2000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human FOXD4L2
Antigen Species:	Human
Gene ID:	286380
Uniprot ID:	Q3SYB3

### Research Background

This gene is a member of the forkhead/winged-helix (FOX) family of transcription factors that have highly conserved FOX DNA-binding domains. Members of the FOX family of transcription factors are regulators of embryogenesis and may play a role in human cancer.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481