

Anti-FOXF1 Antibody (7J721)

Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 40 kDa. Actual: 40 kDa.
Clone:	7J721
Purification:	Protein A purified

Applications

Verified Activity:	4% Paraformaldehyde-fixed Hela (H) cell; Triton X-100 at r. T. for 20 min; Antibody incubation with (FOXF1) monoclonal Antibody, unconjugated (TMAB-06122) 1:100, 90 min at 37°C; followed by conjugated Goat Anti-Rabbit IgG antibody (green) at 37°C for 90 min, DAPI (blue) was used to stain the cell nuclei. PBS instead of the primary antibody was used as the blank control.
Application:	ICC/IF,FCM
Recommended	ICC/IF: 1:50-200; FCM: 1:50-100

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human FOXF1
Antigen Species:	Human
Gene ID:	2294
Uniprot ID:	Q12946

Research Background

Probable transcription activator for a number of lung-specific genes.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481