

## Anti-FRA10AC1 Polyclonal Antibody

## Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Mouse,Rat,Chicken,Dog,Pig,Cow,Horse,Rabbit,Sheep)
Molecular Weight:	Theoretical: 38 kDa. Actual: 38 kDa.
Purification:	Protein A purified

## Applications

Verified Activity:	<p>1. Paraformaldehyde-fixed, paraffin embedded (Human kidney); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30 min; Antibody incubation with (C10orf4) Polyclonal Antibody, Unconjugated (TMAB-06155) at 1:400 overnight at 4°C, followed by operating according to SP Kit (Rabbit) instructions and DAB staining.</p> <p>2. Sample: 293T (Human) Cell Lysate at 40 µg Primary: Anti-C10orf4 (TMAB-06155) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 38 kD Observed band size: 38 kD</p>
Application:	WB,IHC-P,IHC-Fr,IF
Recommended	WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human C10orf4
Antigen Species:	Human
Gene ID:	118924

## Research Background

The protein encoded by this gene is a nuclear phosphoprotein of unknown function. This gene contains a tandem CGG repeat region within a CpG island that normally consists of 8-14 repeats but can expand to over 200 repeats. The repeat region is within the 5' UTR of some transcript variants, but is intronic to another variant. The expanded repeat allele is a fragile site and becomes hypermethylated, causing a reduction in gene expression. A disease phenotype has not been associated with expanded alleles. This gene is found within the rare FRA10A folate-sensitive fragile site. [provided by RefSeq, Dec 2016]

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481