

Anti-FTMT Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Mouse,Rat,Pig,Cow,Horse,Rabbit)
Molecular Weight:	Theoretical: 22 kDa. Actual: 25 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	Western blot analysis of lysates from VEC cells, primary antibody was diluted at 1:1000, 4°over night
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human FTMT
Antigen Species:	Human
Gene ID:	94033
Uniprot ID:	Q8N4E7

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481