

Anti-GABPA Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Rat (predicted:Mouse,Chicken,Dog,Pig,Cow,Rabbit,Sheep)
Molecular Weight:	Theoretical: 51 kDa. Actual: 51 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	<p>1. Paraformaldehyde-fixed, paraffin embedded (Rat kidney); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30 min; Antibody incubation with (GABPA) Polyclonal Antibody, Unconjugated (TMAB-06252) at 1:400 overnight at 4°C, followed by operating according to SP Kit (Rabbit) instructions and DAB staining.</p> <p>2. Sample: 293T Cell (Human) Lysate at 30 µg Primary: Anti-GABPA (TMAB-06252) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 51 kD Observed band size: 51 kD</p>
Application:	WB,IHC-P,IHC-Fr,IF
Recommended	WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human GABPA
Antigen Species:	Human
Gene ID:	2551
Uniprot ID:	Q06546

Research Background

The transcription factor GA-binding protein (GABP) is composed of two subunits, the Ets-related GABP-alpha and a GABP-alpha-associated subunit, GABP beta. GABP alpha binds to a specific DNA sequence and GABP beta exists as b1 and b2 splice variants that differ in their C-termini. In primary neuronal cultures, GABP beta is expressed in both the cytoplasm and the nucleus, whereas GABP alpha is expressed mainly in the nucleus. GABP is constitutively expressed as either a GABP alpha beta heterodimer or a GABP alpha b heterotetramer, both of which can modify GABP-dependent transcription in vitro and in vivo. The GABP alpha beta tetrameric complex performs many different functions, such as stimulating transcription of the adenovirus E4 gene, differentially activating BRCA1 expression in human breast cell lines, potentiating Tat-mediated activation of long terminal repeat promoter transcription and viral replication in certain cell types, acting as a coordinator of mitochondrial and nuclear transcription for

cytochrome oxidase in neurons and assisting in the regulation of rpL32 gene transcription.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481