

Anti-Galactosidase alpha Antibody (1Z79)

Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 49 kDa. Actual: 47 kDa.
Clone:	1Z79
Purification:	Protein A purified

Applications

Verified Activity:	<ol style="list-style-type: none">25 µg total protein per lane of various lysates (see on figure) probed with Galactosidase alpha monoclonal antibody, unconjugated (TMAB-06280) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.Paraformaldehyde-fixed, paraffin embedded Human Kidney; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with Galactosidase alpha Monoclonal Antibody, Unconjugated (TMAB-06280) at 1: 200 overnight at 4°C, followed by conjugation to the Goat Anti-Rabbit IgG H&L Secondary Antibody-HRP and DAB staining.Paraformaldehyde-fixed, paraffin embedded Human Adrenal Gland; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with Galactosidase alpha Monoclonal Antibody, Unconjugated (TMAB-06280) at 1: 200 overnight at 4°C, followed by conjugation to the Goat Anti-Rabbit IgG H&L Secondary Antibody-HRP and DAB staining.
Application:	WB,IHC-P,IHC-Fr,ICC/IF,IF,IP,FCM
Recommended	WB: 1:500-2000; IHC-P: 1:50-200; IHC-Fr: 1:50-200; ICC/IF: 1:50-200; IF: 1:50-200; IP: 1:20-50; FCM: 1:50-100

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human alpha Galactosidase A
Antigen Species:	Human
Gene ID:	2717
Uniprot ID:	P06280

Research Background

Catalyzes the hydrolysis of glycosphingolipids and participates in their degradation in the lysosome.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481