

Anti-GALNTL4 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Mouse (predicted:Human,Rat,Chicken,Dog,Cow,Horse,Rabbit,Sheep)
Molecular Weight:	Theoretical: 70 kDa. Actual: 70 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	Sample: Lung (Mouse) Lysate at 40 µg Primary: Anti-GALNTL4 (TMAB-06311) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 70 kD Observed band size: 70 kD
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human GALNTL4
Antigen Species:	Human
Gene ID:	374378
Uniprot ID:	Q6P9A2

Research Background

The UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetylgalactosaminyltransferase (GalNAc-T) family of enzymes are substrate-specific proteins that catalyze the transfer of GalNAc (N-acetylgalactosaminyl) to serine and threonine residues onto various proteins, thereby initiating mucin-type O-linked glycosylation in Golgi apparatus. GalNAc-TL4, also known as LGALS14 or polypeptide GalNAc transferase-like protein 4, is a 607 amino acid single-pass type II membrane protein belonging to the glycosyltransferase 2 family and GalNAc-T subfamily. Localizing to Golgi apparatus, GalNAc-TL4 utilizes manganese and calcium as cofactors and may assist with the transfer of N-acetyl-D-galactosamine to a serine or threonine residue on protein receptors. GalNAc-TL4 likely catalyzes the initial reaction in O-linked oligosaccharide biosynthesis and contains a ricin B-type lectin domain, which binds to GalNAc and contributes to glycopeptide specificity, and two conserved domains located in the glycosyltransferase region. The N-terminal domain, also known as domain A or GT1 motif, may be involved in manganese coordination and substrate binding while the C-terminal domain, also known as domain B or Gal/GalNAc-T motif, is likely involved in catalytic reaction and UDP-Gal binding. GalNAc-TL4 exists as two alternatively spliced isoforms.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481