

Anti-GM130 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Mouse,Rat)
Molecular Weight:	Theoretical: 113 kDa. Actual: 120 kDa.
Purification:	Protein A purified

Applications

	Sample:
	Lane 1: Human MCF-7 cell lysates
	Lane 2: Human HeLa cell lysates
Verified Activity:	Primary: Anti-GM130 (TMAB-06591) at 1/1000 dilution
	Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution
	Predicted band size: 113 kDa
	Observed band size: 120 kDa
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: mouse GM130
Antigen Species:	Mouse
Gene ID:	2801

Research Background

The docking of transport vesicles to their target membrane is mediated by p115. GM130, a cis-Golgi matrix protein, interacts specifically with p115 and provides a membrane docking site. Both GM130 and p115 are involved in vesicle tethering to Golgi membranes. The amino-terminus of GM130 binds to p115, whereas the carboxy-terminus binds to Golgi membranes. Both Giantin and GM130 compete for binding to p115. Thus, p115-Giantin and p115-GM130 interactions might mediate independent membrane tethering events. Transport from the ER to the cis/medial Golgi compartments requires the action of p115, GM130 and Giantin via a sequential rather than a coordinate mechanism. Mitotic phosphorylation of GM130 at Serine 25 is mediated by Cdc2, prevents binding to p115 and is directly involved in mitotic Golgi fragmentation. GM130 is phosphorylated in prophase as the Golgi complex starts to break down, and remains phosphorylated in metaphase and anaphase. In telophase, GM130 is dephosphorylated by PP2A as the Golgi fragments start to reassemble.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481