

Anti-GPATCH1/ECGP Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Mouse,Rat,Dog,Pig,Cow,Horse,Sheep)
Molecular Weight:	Theoretical: 103 kDa. Actual: 103 kDa.
Purification:	Protein A purified

Applications

Sample:	HL60 (Human) Cell Lysate at 30 µg Molt-4 (Human) Cell Lysate at 30 µg
Verified Activity:	Primary: Anti-GPATCH1'ECGP (TMAB-06653) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 103 kD Observed band size: 103 kD
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human GPATCH1/ECGP
Antigen Species:	Human
Gene ID:	55094
Uniprot ID:	Q9BRR8

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481