

Anti-GPR40 Polyclonal Antibody 2

Product Details

Ig Type:	IgG
Reactivity:	Mouse,Rat
Molecular Weight:	Theoretical: 32 kDa. Actual: 28 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	25 µg total protein per lane of various lysates (see on figure) probed with GPR40 polyclonal antibody, unconjugated (TMAB-06715) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: mouse GPR40
Antigen Species:	Mouse
Gene ID:	233081
Uniprot ID:	Q76JU9

Research Background

G protein coupled receptors provide attractive targets for drug therapy due to the sheer size and diversity of ligands within this receptor family. G protein-coupled receptor 40 (GPR40) functions as a cell-surface receptor for long-chain free fatty acids (FFAs). FFAs provide an important energy source, but also function as signaling molecules in various pathways, notably the process of insulin secretion. In pancreatic tissue, the interaction of long chain FFAs with GPR40 amplifies glucose-stimulated insulin secretion from beta cells, suggesting a possible role for GPR40 in the treatment of diabetes associated with insulin-deficiency. Specifically, the Arg211His polymorphism in the GPR40 gene may contribute to the variation of insulin secretory capacity in Japanese men. Also, GPR40 may be involved in the control of breast cancer cell growth by fatty acids and, therefore, provide a link between fat and cancer.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481