

Anti-GSTK1 Antibody (1G279)

Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 25 kDa. Actual: 25 kDa.
Clone:	1G279
Purification:	Protein A purified

Applications

Verified Activity:	Western blot analysis of 293T cell lysate. Using GSTK1 (TMAB-06826) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	IHC-P,IHC-F,ICC/IF,IF,IP,FCM
Recommended	IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:50-200,IF=1:100-500,IP=1:20-50,FCM=1:50-100

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human GSTK1
Antigen Species:	Human
Gene ID:	373156
Uniprot ID:	Q9Y2Q3

Research Background

Glutathione S-transferase that catalyzes the conjugation of glutathione to exogenous and endogenous compounds. Significant glutathione conjugating activity is found only with the model substrate, 1-chloro-2,4-dinitrobenzene.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481