

Anti-HLC8 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Mouse,Rat,Dog,Cow,Sheep)
Molecular Weight:	Theoretical: 67 kDa. Actual: 67 kDa.
Purification:	Protein A purified

Applications

Sample:	Jurkat (Human) Cell Lysate at 30 µg K562 (Human) Cell Lysate at 30 µg HL60 (Human) Cell Lysate at 30 µg
Verified Activity:	Primary: Anti-HLC8 (TMAB-07172) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 67 kD Observed band size: 67 kD
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human HLC8/Lung cancer related protein 8
Antigen Species:	Human

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481