

Anti-HLTF Antibody (1B226)

Product Details

| | |
|-------------------|--|
| Ig Type: | IgG |
| Reactivity: | Human |
| Molecular Weight: | Theoretical: 114 kDa. Actual: 130 kDa. |
| Clone: | 1B226 |
| Purification: | Protein A purified |

Applications

| | |
|--------------------|--|
| Verified Activity: | <p>1. 25 ug total protein per lane of various lysates (see on figure) probed with HLTF monoclonal antibody, unconjugated (TMAB-07173) at 1: 1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.</p> <p>2. The HeLa (H) cells were fixed with 4% PFA (10 min at r. T.) and then permeabilized with 90% ice-cold methanol for 20 min at -20°C, the cells then were incubated in 5% BSA to block non-specific protein-protein interactions (30 min at r. T.), followed by secondary antibody incubation for 40 min at room temperature. Primary Antibody (green): Rabbit Anti-HLTF antibody (TMAB-07173, 1: 100); Isotype Control (orange): Rabbit IgG. Blank control (black): PBS. Acquisition of 20,000 events was performed.</p> |
| Application: | WB,FCM |
| Recommended | WB=1:500-2000,FCM=1:50-100 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|-----------------------------------|
| Immunogen: | A synthesized peptide: human HLTF |
| Antigen Species: | Human |
| Gene ID: | 6596 |
| Uniprot ID: | Q14527 |

Research Background

Has both helicase and E3 ubiquitin ligase activities. Possesses intrinsic ATP-dependent nucleosome-remodeling activity; This activity may be required for transcriptional activation or repression of specific target promoters (By similarity).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481