

Anti-human CD86 Antibody (5D723)

Product Details

| | |
|-------------------|---|
| Ig Type: | Mouse IgG1, k |
| Reactivity: | Human |
| Molecular Weight: | Theoretical: 31 kDa. Actual: 72-74 kDa. |
| Clone: | 5D723 |
| Purification: | Protein G purified |

Applications

Blank control: THP-1.
Primary Antibody (green line): Mouse Anti-CD86 antibody (TMAB-07422)
Dilution: 10 ul/Test;
Secondary Antibody (white blue line) : Goat anti-Mouse IgG-FITC
Dilution: 0.5 µg/Test.

Verified Activity: Isotype control (orange line): Normal Mouse IgG
Protocol

The cells were incubated in 5% BSA to block non-specific protein-protein interactions for 30 min at room temperature. Cells stained with Primary Antibody for 30 min at room temperature. The secondary antibody used for 40 min at room temperature. Acquisition of 20,000 events was performed.

Application: FCM

Recommended FCM: 10µl/Test

Properties

Stability & Storage: Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Gene ID: 942

Research Background

This gene encodes a type I membrane protein that is a member of the immunoglobulin superfamily. This protein is expressed by antigen-presenting cells, and it is the ligand for two proteins at the cell surface of T cells, CD28 antigen and cytotoxic T-lymphocyte-associated protein 4. Binding of this protein with CD28 antigen is a costimulatory signal for activation of the T-cell. Binding of this protein with cytotoxic T-lymphocyte-associated protein 4 negatively regulates T-cell activation and diminishes the immune response. Alternative splicing results in several transcript variants encoding different isoforms.[provided by RefSeq, May 2011].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481