

Anti-JSRV-ENV Antibody (2Q457)

Product Details

| | |
|-------------------|--------------------------------------|
| Ig Type: | IgG |
| Reactivity: | JSRV |
| Molecular Weight: | Theoretical: 66 kDa. Actual: 37 kDa. |
| Clone: | 2Q457 |
| Purification: | Protein G purified |

Applications

| | |
|--------------------|---|
| Sample: | Lane 1: Recombinant JSRV-ENV protein, His |
| Verified Activity: | Primary: Anti-JSRV-ENV (TMAB-07885) at 1/1000 dilution |
| | Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution |
| | Predicted band size: 66 kDa Observed band size: 37 kDa |
| Application: | WB |
| Recommended | WB: 1:500-2000 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------|-------------------------------|
| Immunogen: | Recombinant Protein: JSRV-ENV |
|------------|-------------------------------|

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481