

Anti-JSRV-ENV Antibody (2Q457)

Product Details

Ig Type:	IgG
Reactivity:	JSRV
Molecular Weight:	Theoretical: 66 kDa. Actual: 37 kDa.
Clone:	2Q457
Purification:	Protein G purified

Applications

Verified Activity:	Sample:
	Lane 1: Recombinant JSRV-ENV protein, His
	Primary: Anti-JSRV-ENV (TMAB-07885) at 1/1000 dilution
	Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution
	Predicted band size: 66 kDa
	Observed band size: 37 kDa
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: JSRV-ENV
------------	-------------------------------

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481