

Anti-KBTBD2 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse (predicted:Rat,Dog,Pig,Cow,Horse)
Molecular Weight:	Theoretical: 71 kDa. Actual: 71 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	1. Sample: U251 (Human) Cell Lysate at 30 µg Primary: Anti-KBTBD2 (TMAB-07915) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 71 kD Observed band size: 71 kD
	2. Sample: Cerebrum (Mouse) Lysate at 40 µg Primary: Anti-KBTBD2 (TMAB-07915) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 71 kD Observed band size: 71 kD
	Application: WB
	Recommended WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human KBTBD2
Antigen Species:	Human
Gene ID:	25948
Uniprot ID:	Q8IY47

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481