

Anti-KCNK3 Polyclonal Antibody 2

Product Details

Ig Type:	IgG
Reactivity:	Mouse,Rat (predicted:Human,Chicken,Dog,Cow,Rabbit,Sheep)
Molecular Weight:	Theoretical: 43 kDa. Actual: 43 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	1. Sample: Pancreas (Mouse) Lysate at 40 µg
	Primary: Anti-KCNK3 (TMAB-07940) at 1/300 dilution
	Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution
	Predicted band size: 43 kD
	Observed band size: 43 kD
	2. Sample: Placenta (Mouse) Lysate at 40 µg
	Primary: Anti-KCNK3 (TMAB-07940) at 1/300 dilution
	Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution
	Predicted band size: 43 kD
Observed band size: 43 kD	3. Sample: Urinary bladder (Mouse) Lysate at 40 µg
	Primary: Anti-KCNK3 (TMAB-07940) at 1/1000 dilution
	Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution
Predicted band size: 43 kD	
Observed band size: 49 kD	
Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human KCNK3
Antigen Species:	Human
Gene ID:	3777
Uniprot ID:	O14649

Research Background

KCNK9 or TASK-3 (TWIK-related Acid sensitive K⁺ channel) is a member of the potassium channel family of proteins that contain two-pore domain and four transmembrane domains. These channels are characterized as leak K⁺ channels that are sensitive to changes in the extracellular pH. The physiological functions of TASK channels are largely unknown; it has been proposed that they may be involved in the regulation of breathing, aldosterone

A DRUG SCREENING EXPERT

secretion and anesthetic-mediated neuronal activity. They were found to act in neurons' membrane potential and in resting K⁺ currents. KCNK9 gene has been localized to the chromosomal region 8q24. The protein is expressed at high levels mainly in the brain and at low levels in other tissues. In contrast to normal tissues, it was found that KCNK9 is amplified in some human carcinomas such as breast, lung, colon and metastatic prostate. In 10% of breast cancer patients this gene is amplified, and in 44% the protein is over expressed. Monoclonal antibodies to KNCK9 are an important tool for studying the potassium channel family of proteins in different tissues.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481