

Anti-KIAA0649 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human (predicted:Mouse,Rat,Dog,Horse)
Molecular Weight:	Theoretical: 127 kDa. Actual: 127 kDa.
Purification:	Protein A purified

Applications

Verified Activity:	1. Sample: A549 (Human) Cell Lysate at 40 µg Primary: Anti-KIAA0649 (TMAB-08004) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 127 kD Observed band size: 127 kD
	2. Sample: DU145 (Human) Cell Lysate at 40 µg Primary: Anti-KIAA0649 (TMAB-08004) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 127 kD Observed band size: 127 kD
	Application: WB
	Recommended WB: 1:500-2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human KIAA0649
Antigen Species:	Human
Gene ID:	9858
Uniprot ID:	Q5T8A7

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481