

## Anti-MAS1L Antibody (5W447)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Molecular Weight:	Theoretical: 42 kDa. Actual: 37 kDa.
Clone:	5W447
Purification:	Protein A purified

### Applications

Verified Activity:	Western blot analysis of 293T cell lysate. Using MAS1L (TMAB-08623) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r. T. for 60 min.
Application:	WB,IP
Recommended	WB: 1:1000-2000; IP: 1:20-50

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human MAS1L
Antigen Species:	Human
Gene ID:	116511
Uniprot ID:	P35410

### Research Background

Belongs to the G-protein coupled receptor 1 family. Mas subfamily.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481